

**Tool ID** 15J1359AA/15J1358AA/15J1357AA

**a. General**

Cat	Description
<b>Vendor / Supplier</b>	Waters 2695 Separation Module and other components
<b>Model</b>	2695
<b>Process</b>	GPC - Gel Permeation Chromatography
<b>Vintage (Year of Manufacture)</b>	Unknown
<b>Serial No.</b>	J00SM4
<b>Asset Description</b>	Gel Permeation Chromatography system with Photodiode array detector and refractive index detector
<b>Tool Status</b>	Cold shutdown
<b>Software Version</b>	MS Windows xp
<b>CIM</b>	x

**b. Hardware Configuration (Fab)**

Cat	Description	Qty / Units	Status,tick if ok[√ ]
<b>Components</b>	Waters 2695	1	ok[√ ]
	Waters 2410	1	ok[√ ]
	Water 996	1	ok[√ ]
	Computer	1	ok[√ ]
	Temperature Control Module	1	ok[√ ]

**c. Hardware Configuration (Subfab / Auxilliary Units)**

Sub Accessories / Auxilliary Units	Qty / Units	Status,tick if ok[√ ]
NA		

**d. Missing/Faulty Parts / Accesories List**

**(If NO missing/faulty parts/accessories please indicate NA)**

Description	Qty / Units
NIL	

Tool ID : 15J1359AA/15J1358AA/15J1357AA.

**a. General**

Description	Photos
Eqpt ID	
Vendor / Supplier / Model	15J1359AA/15J1358AA/15J1357AA

**b. Hardware Configuration (Fab)**

Description	Photos
Main System	  
	